



Drug of Abuse



Fast Result
in 5 minutes

99.86%
Accuracy

99.93%
Specificity

99.80%
Sensitivity





Easy to Use, Easy to Read Convenient Visual Interpretation

Drug abuse is a disorder that is characterized by a destructive pattern of using a substance that leads to significant problems or distress.

Teens are increasingly engaging in prescription drug abuse, particularly narcotics, and stimulant medications, which treat conditions like attention deficit disorder and narcolepsy.

TEST PROCEDURE

1. Device



5 Minutes



Positive




Negative



Invalid

2. Panel



5 Minutes



Positive



Negative



Invalid

3. Strip



5 Minutes



Positive



Negative



Invalid

PARAMETER	SPECIMEN	ACCURACY (ng/ml)	ISI/BOX STRIP / CARD
1. Amphetamine	Urine	1000 ng/ml	50 / 25
2. Methamphetamine	Urine	300 ng/ml	50 / 25
3. Benzodiazepine	Urine	300 ng/ml	50 / 25
4. Cocaine	Urine	1000 ng/ml	50 / 25
5. Marijuana (THC)	Urine	50 ng/ml	50 / 25
6. Opiate	Urine	300 ng/ml	50 / 25
7. Barbiturate	Urine	300 ng/ml	50 / 25
8. Ketamine	Urine	1000 ng/ml	50 / 25
9. MDMA	Urine	500 ng/ml	50 / 25
10. Methadone	Urine	300 ng/ml	50 / 25

MULTIPLE DRUG PANEL TEST	SPECIMEN	BOX
1. DOA 3 (AMP/MOR/THC)	Urine	25 Cards
2. DOA 5 (AMP/MOR/THC/BZD/MET)	Urine	25 Cards
3. DOA 6 (AMP/MOR/THC/BZD/MET/COC)	Urine	25 Cards

